

TestRoom

VISUAL TRAINING

PROFESSIONAL ARTICLES
FOR EYESIGHT TESTING



7.0 REFRACTION & ORTHOPTIC

Our products comply with **Regulation (EU) no. 2017/745**.

Regulation (EU) no. 2017/745 on medical devices establishes the general safety and performance requirements of devices with the aim of guaranteeing the safety and protection of the health of patients, users and any third parties during the use of such products.

MD MD Medical device compliant with EU Regulation no. 2017/745



7.1 Home Visits

- Home Visits 4
- Suitcase for eye check / Home care 5

7.2 Press-On - Occluders

- Strabismus and Amblyopia 6
- Press-On Prism 8
- Press-On Lenses 10
- Bangerter Occlusion Foils 12
- Hypoallergenic occluders for frames 14
- Hypoallergenic occluders with elastic 15
- Adhesive Eye Patches 16

7.3 Visual Training

- Introduction 18
- Marsden Spheres 20
- Near Point Dot Card - Brock String 21
- Red/Green reading sheet and vinyl sheets 22
- Hendrickson lifesaver card 23
- Red & Green wraparound glasses 23
- Pairs of cards, transparent and opaque, with eccentric red/green circles and polarized 24
- Binocular vision assessment test 25
- Kit aperture rule 26
- Bull's Eye - Hart's Tables 27
- Wolff's double sphere target 28
- Lang Fixation Cubes - Lang Fixation Sticks 29
- Flip-Up occluder glasses 30
- Glasses with R/G filter - Tiger shape 31
- Glasses with polarized lenses - Parrot shape ... 31
- Home Trainer Kit 32
- Lectern 34



It will show you the items that can be used by users at home to carry out visual training exercises, after an optometric check.



7.4 Diagnostics

• The beginning of visual analysis	36
• Berens Prism Bar	38
• Retinoscopy Racks	39
• Test for determining eye fixation dominance ..	40
• Ishihara Colour Test	41
• Prism 6Δ	42
• White Light Pen - Worth 4 dot flashlight	43
• Hand held occluders	44
• Fly Stereopsis Test - Butterfly Stereopsis Test ..	46
• Random Dot Stereopsis Test - Polarized Goggles ..	47
• Lang Stereo Test® I-R	48
• Teddy Bear Test	49
• Amsler Test	50

7.5 Flippers: cross cylinders and confirmation test

• Metal flipper binocular confirmation test	51
• Metal flipper cross cylinders	51
• Plastic flippers binocular confirmation test	52
• Plastic flippers polarized and red-green test ...	52
• Plastic prismatic flippers	53

7.6 Trial Frames - Trial Lens Set

• Children Trial Frame Kit	54
• Trial Frames	56
• Meniscus Trial Lens Set	62
• Regular Trial Lens Set	63
• Progressive Trial Lens Set	66
• Prism Trial Lens Set	67
• Trial Clips and Trial Lens Holder	68

7.7 Equipment and accessories

• Digital PD Meter	69
• Refractor/Phoropter - Vision Tester	70
• Self-balancing Arm for Refractor/Phoropter ...	71
• Multifunction Ruler	72
• PD Scales - Adhesive film	73
• Reading charts distance tests	74
• Optotypes for near vision	76

7.8 Contact Lens Accessories - Poster

• Schirmer Test	77
• Fluorescein Sodium	77
• Easy Mark	77
• Cases For soft and G.P. hard lenses	78
• Suckers Holders and Tweezers	79
• Disposable tissues	80
• Laminated poster "THE EYE"	81

7.1 HOME VISITS

Service, relationship, business opportunity

Centro Style offers a range of useful products suitable for visits outside the optical centre. An increasingly requested and appreciated service, an important moment of relationship with customers who, due to problems relating, for example, to advancing age and travel, do not have the possibility of going to the store to have their eyesight checked. Easy portability, simplicity of use of the instruments (both classic and innovative) and minimal bulk are some of the prerogatives of this system, which aims to make the home visit effective and complete.



- *Service*
- *Innovation*
- *Professionalism*
- *Relationship*
- *Differentiation*
- *Opportunities*

Suitcase for eye check/home care

Designed to hold a selection of items ideal for an in-home visual exam or for screenings outside the office. It also contains eyewear maintenance and adjustment tools consisting of a selection of screws, nose pads, pliers and screwdrivers most used in daily practice.



TABLET NOT INCLUDED

NEW

Ref. 08269 **Kr 13 700,-**

Contains:

For eyesight testing

- Ref. 08262 Regular trial lens set with 90 lenses
- Ref. 08311 Red occluder
- Ref. 08312 Frosted occluder
- Ref. 08313 Black occlude
- Ref. 08317 Pen light
- Ref. 08333 Confirmation Test Visus +/-0.50D
- Ref. 04961 Multifunction ruler
- Ref. 04959 Plastic PD ruler with black occluder
- Ref. 08427 Reading chart 3 m distance with "E" letters

Dimensions: 555x370x95 mm

Weight: 5,5 kg.

Tablet not included

For the maintenance and adjustment of the glasses:

- Ref. 03250N Easy Move inclination plier
- Ref. 03722 Screw holding tweezer
- Ref. 02257 Ball-bearing screwdriver
- Ref. 03546N Hard metal inclined side cutter
- Ref. 03200N Easy move holding plier
- Ref. 00299 Stainless steel screws for flex hinges
- Ref. 00410 Self-tapping stainless steel screws M1.2
- Ref. 00412 Self-tapping stainless steel screws M1.4
- Ref. 00341 Nickel-silver screws for nose pads M1,0
- Ref. 01466 Symmetrical screw fixing silicone nose pads 14,5 mm
- Ref. 01461 Symmetrical push-on silicone nose pads 14,5 mm

7.2 VISUAL TRAINING

Strabismus and Amblyopia

Specific pathologies or visual problems, if identified early (within six years of age), can be effectively improved until completely resolved thanks to specific and timely visual rehabilitation interventions.

Among the many ocular anomalies that can develop in the visual system in childhood, the most frequent and important are strabismus and amblyopia.

STRABISMUS

Strabismus is a vision condition in which the eyes are not aligned correctly and point in different directions. In a normal state of vision, both eyes look in the same direction, allowing the brain to receive two separate images and combine them into one three-dimensional image.

However, in people with strabismus, one of the eyes may turn inward (convergent strabismus), outward (divergent strabismus), up or down, causing a mismatch

between the points of the object and points in the eye that can cause double vision (diplopia) or suppression of the image from the deviated eye, to avoid double vision.

Strabismus can be present from birth or develop in childhood. If left untreated, it can lead to vision problems and development of three-dimensional perception.

Alterations in binocular vision induce a series of sensory adaptations, aimed at minimizing the disadvantages of the deviation. Three types of adaptation are identified based on the age at which the deviation appears and its severity:

- monofixation syndrome if the angle of deviation is very small, <10DP;
- anomalous retinal correspondence if the angle of deviation is between 10 and 20DP;
- large suppression scotoma if the angle of deviation is very large, greater than 20DP (Nucci, Serafino 2012).



The rehabilitation path to undertake depends on the causes that generate diplopia (double vision). Regardless of the cause, it is possible to address the problem with techniques that allow you to eliminate visual duplication. In particular it is possible to use special adhesive lenses that have prismatic power called Press-On prism or additional power called Press-On lens.

Strabismus and Amblyopia

AMBLYOPIA

Amblyopia is a condition characterised by a reduction in the visual capacity of one or, more rarely, both eyes, which can be caused by various eye diseases. The word 'amblyopia' comes from the Greek and means 'dim vision', thus defining the deficit in the maturation process of the visual pathways induced by abnormal visual stimulation.

Clinically, an eye is defined as amblyopic if it has a visual acuity that is at least two-tenths less than the other eye; furthermore, the amblyopic eye also presents a reduction in contrast sensitivity and altered chromatic perception (Rajavi et al., 2015; Kocakaltintas et al., 2000). The amblyopic eye is often called lazy because it is not normally used by the brain and therefore has a visual capacity that it does not actually use.

Amblyopia can occur in strabismus, in the presence of very high or significantly different refractive errors

between the two eyes (anisometropia) and finally in the presence of obstacles on the visual axis (congenital cataract, eyelid ptosis, etc.) which prevent the light stimulus to arrive on the retina in an adequate manner. Amblyopia can also arise in the presence of micro strabismus, i.e. in the case in which the angle of deviation is very small; in this case, double vision of objects does not occur.

This visual defect, if treated early and adequately in preschool age, is largely recoverable, otherwise it becomes a visual condition that will persist even into adulthood when it can no longer be modified.

Rehabilitation of amblyopia is carried out by **occluding** or **penalizing** the eye with the **best vision** - by means of plasters or occluders to be applied directly to the eye or fabric occluders or suction cups with or without elastic to be applied to the glasses - **to stimulate the lazy eye** to take action to improve their visual ability.



7.2 PRESS-ON PRISM

For the treatment of strabismus

Kr 478,-



Prism 1^Δ



Prism 10^Δ



Prism 40^Δ

PRESS-ON™ PRISM: Used to treat several ocular motility disorders including strabismus.

Advantage of Press-On: practical, lightweight, wide power range. Thin and flexible and easy to fix to existing lenses - held by static adhesion, easily removed. Low cost and low risk treatment during power changes and vision therapy .

PRISM POWERS IN DIOPTRES - Prism 67 mm - 1 pc.

Ref.	Power	Ref.	Power	Ref.	Power
08101	1 ^Δ	08107	7 ^Δ	08116	20 ^Δ
08102	2 ^Δ	08108	8 ^Δ	08117	25 ^Δ
08103	3 ^Δ	08109	9 ^Δ	08118	30 ^Δ
08104	4 ^Δ	08110	10 ^Δ	08119	35 ^Δ
08105	5 ^Δ	08112	12 ^Δ	08120	40 ^Δ
08106	6 ^Δ	08115	15 ^Δ		

The pictures show how the image is moved by the prism to accommodate the different direction of the patient's line of sight. The image is moved in the opposite direction to the base.

HOW TO APPLY PRESS-ON PRISMS

- There is no need to determine the optical centre of the lens.
- The glossy side of the Press-On is applied to the inside of the lens.



1. Place the lens onto the Press-On with the inner side of the lens facing the smooth surface of the prism, according to the prescription.
2. Trace the shape of the lens onto the Press-On.
3. Cut the Press-On (1 mm less than the traced shape).
4. Moisten the inside of the lens and apply the Press-On prism with the glossy side facing the inside surface of the lens. Do not use any adhesives.



CORRECT POSITIONING OF PRESS-ON PRISM ON LENS

The base orientation of the prism is indicated with the word "base" on the prism edge. Mark lens and prism at the same position on their edges to indicate alignment. The Press-On must be positioned onto the lens according to the BASE description indicated in the prescription: NASAL, TEMPORAL, TOP, BOTTOM, DIAGONAL.

VIDEO



Press-On and 3M are a registered trade mark of 3M

MD MD Medical device compliant with EU Regulation no. 2017/745

For the treatment of strabismus




Vertical Correction

NASAL BASE

The prism's lines are vertical.
 Real position of the word 'BASE' on the Press-On:
 Right eye: nasal, at the bottom, text vertical.
 Left eye: nasal, at the top, text vertical.

TEMPORAL BASE

The prism's lines are vertical.
 Real position of the word 'BASE' on the Press-On:
 Right eye: temporal, at the top, text vertical.
 Left eye: temporal, at the bottom, text vertical.



Horizontal Correction


TOP BASE

The prism's lines are horizontal.
 Real position of the word 'BASE' on the Press-On:
 Right eye: nasal, at the top, text horizontal.
 Left eye: temporal, at the top, text horizontal.

BOTTOM BASE

The prism's lines are horizontal.
 Real position of the word 'BASE' on the Press-On:
 Right eye: temporal, at the bottom, text horizontal.
 Left eye: nasal, at the bottom, text horizontal.

Oblique Correction



DIAGONAL BASE

The prism's lines are diagonal.
 Such a prescription is achieved by locating the axis of the prism diagonally determining the angle with a Rotation Nomograph.



7.2 PRESS-ON

Lenses



Press-On and 3M are a registered trade mark of 3M



CE Regulation (EU) 2017/745



Kr 478,-

PRESS-ON™ LENS :

Advantage of Press-On: practical, lightweight, widerange of powers. Thin and flexible and easy to fix to existing lenses - held by static adhesion, easily removed. Low cost and low risk treatment during power changes and vision therapy.

POSITIVE LENSES - Lens 67 mm - 1 pc.

Ref.	Power	Ref.	Power
08121	+1	08123	+2,5
08162	+1,5	08124	+3
08122	+2	08126	+4

HOW TO APPLY PRESS-ON LENSES

1. Determine the optical centre of the lens and "dot" it on the front surface.
2. Place the lens onto the Press-On with the inner side of the lens facing the smooth surface of the Press-On.
3. Make sure that the optical centre of the lens coincides with the optical centre of the Press-On, which is the smallest concentric ring.
4. Trace the shape of the lens onto the Press-On.
5. Cut the Press-On, but slightly smaller than the lens.
6. Moisten the inside of the lens and apply the Press-On with the glossy side facing the inside surface of the lens.

VIDEO



MD MD Medical device compliant with EU Regulation no. 2017/745

KIT PRESS-ON

7.2


To apply Press-On prism and Press-On Lenses

Ref. 08139 **Kr 729,-**

Kit consisting of:

- Fine tip marker (Ref. 06970)
- Scissor for shaping (Ref. 03782)
- Protractor to apply Press-On with oblique base (Ref. 08135)
- Tray (Ref. 01667)

VIDEO



PROTRACTOR

Ref. 08135 **Kr 132,-**


with a correct inclination. 2 pcs.



SCISSOR

Ref. 03782 **Kr 189,-**

To cut and shape Press-On.



MARKER

Ref. 06970 **Kr 232,-**

Ideal to mark the shape of the lens.

7.2 BANGERTER

Graded foils - Occlusion foils

Kr 109,-


Occlusion foils are a system of graded thin flexible foils of varying degrees of transparency that equalize the spatial contrast of the dominant eye to that of the amblyopic eye. They cling to the spectacle lens of the healthy eye to improve the function of the weaker eye. There are 10 levels of occlusion - see table below.

Use: These are used to treat amblyopia, exclusion and the anomalous retinal correspondence found in amblyopia; problems relating to correlation, monolateral aphakia and temporary diplopia.




Examples of vision with different occlusion levels


Dimensions: 65x65 mm. - 1 pc. / Ref.




Ref. 08001 Visus ~1.0 - N 1
Visual Acuity ~20/20




Ref. 08002 Visus ~0.8 - N 2
Visual Acuity ~20/265




Ref. 08003 Visus ~0.6 - N 3
Visual Acuity ~20/30




Ref. 08004 Visus ~0.4 - N 4
Visual Acuity ~20/50




Ref. 08005 Visus ~0.3 - N 5
Visual Acuity ~20/70



Ref. 08009 Visus ~0.2 - N -
Visual Acuity ~20/100




Ref. 08006 Visus ~0.1 - N 6
Visual Acuity ~20/200



Ref. 08007 Visus <math><0.1</math> - N 7
Visual Acuity ~20/300




Ref. 08000 Visus LP - N 0
Visual Acuity none



Ref. 08008 Visus 0.0
Visual Acuity none

Graded foils - Occlusion foils



Ref.	08008	08000	08007	08006	08009	08005	08004	08003	08002	08001
Vision	0.0	Light perception	<0.1	0.1	0.2	0.3	0.4	0.6	0.8	1.0
Corr. No.	8	0	7	6	-	3	4	5	2	1
Approx Visual Acuity	none	none	~20/300	~20/200	~20/100	~20/70	~20/50	~20/30	~20/25	~20/20

Ref. 08017 Kr 1231,-
 Kit of assorted Bangerter occlusion foils with display showing reference and 10 levels of occlusion.
 QR code for video.

Ref. 11315 Kr 136,-
 Display only
 Dimensions: 240x240x10 mm



VIDEO



HOW TO APPLY BANGERTER OCCLUSION FOILS

The shiny side of the foil is applied to the inside of the lens.

- Cut the foil to a shape slightly smaller than that of the lens
- Clean the lens thoroughly
- Damp the shiny side of the foil with warm water and apply to the inner side of the lens.
- Press with a dry cloth to remove any air bubbles and traces of damp.
- Allow to dry for 1 hour.



MD MD Medical device compliant with EU Regulation no. 2017/745

7.2 HYPOALLERGENIC OCCLUDERS FOR FRAMES

For orthoptic therapies of amblyopia and strabismus - In fabric

BABY

Dimensions:
105x50 mm
1 pc.

Ref. 08174
Traffic
Kr 212,-

Ref. 08175
Bunnies
Kr 212,-

Ref. 08176
Teddy
Kr 212,-

Ref. 08177
Ducklings
Kr 212,-



Occluders for children and adults frames made of high quality ecological fabric, reusable and washable.

They are made with certified materials, hypoallergenic and breathable, and eco-friendly. The fabric is printed using inks and water based dyes.

They are designed to comfortably fit on any type of frame and provide a total occlusion, effective for the treatment of amblyopia (lazy eye).

These occluders offer the perfect combination of comfort, style and fit to make treatment time fun and effective for children, and elegant for adults.

These occluders have colorful and fun drawings graphics suitable for both the right and left eye.

CHILDREN

Dimensions:
118x60 mm

Ref. 08180
Soccer Blues
Kr 216,-

Ref. 08181
Pirate Skulls
Kr 216,-

Ref. 08182
Jurassic
Kr 216,-

Ref. 08183
Fairies
Kr 216,-

Ref. 08184
Hearts&Flowers
Kr 216,-

Ref. 08185
Hoots
Kr 216,-

NEW **Ref. 08186**
Fire Engines
Kr 216,-

NEW **Ref. 08187**
Flutterby
Kr 216,-



ADULT

Dimensions: 130x60 mm
1 pc.



Ref. 08171 **Kr 259,-**
Regular Beige



Ref. 08172 **Kr 259,-**
Regular Black

MD MD Medical device compliant
with EU Regulation no. 2017/745

HYPOALLERGENIC OCCLUDERS WITH ELASTIC 7.2

For orthoptic therapies of amblyopia and strabismus - In fabric

Occluders for children and adults with elastic to place directly on the eye, usable even without frames.

Entirely made of breathable certified fabric, safe in contact with the skin to ensure a soft and enveloping fit around the eye. Reusable and washable.

Ideal for occluding one eye: they are kept in position with a quality elastic, latex-free, which makes them comfortable to wear for long periods.

The occluders have a black lining to exclude most of the light and a reinforced fabric core that keeps its shape even after the washing.

These occluders have colorful and fun drawings graphics suitable for both the right and left eye.



ADULT
Dim.: 77x67 mm
Ref. 08170 **Kr 384,-**
Beige 1 pc.

EYE PATCH WITH ELASTIC CORD
Ref. 08315 **Kr 342,-**
Total occlusion, in vinyl
Dimensionis: 76,5x52,5 mm
3 pcs.



CHILDREN
Dim.: 70x55 mm
1 pc.

Ref. 08198
Skulls
Kr 216,-

NEW
Ref. 08194
Mighty Monster
Kr 216,-

NEW
Ref. 08195
T-Rex
Kr 216,-

Ref. 08196
Football Crazy
Kr 216,-

Ref. 08197
Hearts
Kr 216,-

Ref. 08199
Magic Unicorn
Kr 216,-

Ref. 08321 **Kr 337,-**
Translucent, in soft plastic
Dim. 76,5x52,5 mm.
3 pcs.



7.2 ADHESIVE EYE PATCHES

Hypo-allergenic fun patches for children's vision problems - **Latex Free**

Reserved for the Italian Market



Each pack contains funny patch stickers.

Each box contains 20 eye patches sealed in individual sleeves.

JUNIOR WHITE
67x50 mm
Ref. 08076
6 Packs
Kr 902,-

MEDIUM WHITE
76x54 mm
Ref. 08079
6 Packs
Kr 902,-

REGULAR WHITE
85x59 mm
Ref. 08077
6 Packs
Kr 820,-



Each box contains 20 eye patches sealed in individual sleeves.

JUNIOR SYMPATHY
67x50 mm
Ref. 08071
1 Pack
Kr 202,-

MEDIUM SYMPATHY
76x54 mm
Ref. 08072
1 Pack
Kr 202,-

- Made of a perforated monostretch non-woven material for extreme comfort.
- A special black screen in the central pad prevents the passage of light. For this reason, they are suitable to treat strabismus and amblyopia.
- Totally hypo-allergenic, latex free, for sensitive skin around the eye. Well ventilated to avoid sweating.
- Can be used for prolonged treatment.
- Close eye when applying.



MD MD Medical device compliant with EU Regulation no. 2017/745

3M HYPO ALLERGENIC EYE PATCHES

7.2

For orthoptic therapies of amblyopia and strabismus - **Latex Free**



Each pack contains 30 patches.

**OPTICLUDE
BOY&GIRLS MINI**
50x62 mm
Ref. 08069
1 Pack
Kr 228,-

**OPTICLUDE
BOY&GIRLS MAXI**
57x82 mm
Ref. 08073
1 Pack
Kr 228,-

- Made of perforated non-woven material and a pad placed on a transpirable and conformable support, covered by a hypoallergenic adhesive.
- Well ventilated to respect the physiological functions of the skin.
- Totally safe for the skin of the periocular area and can be used for prolonged treatments.



- Make sure that the skin around the eye is dry and clean and that the child's face is relaxed before applying the eye patch over the closed eye.
- Open the wrapper and peel off the film from the adhesive side of the patch.
- Apply the patch to the skin with the smaller corner nearest the nose.
- Use fingers to apply pressure to all outer edges of the eye patch to ensure that all of the patch is sticking to the skin.
- We recommend to replace eye patches every day or according to the doctor's instructions.



MD MD Medical device compliant with EU Regulation no. 2017/745

7.3 VISUAL TRAINING

Exercises to improve the effectiveness of visual skills

Visual training is an Optical Optometric activity aimed at promoting the development, improvement and optimization of the visual process. It is a set of specific exercises and procedures that train and stimulate visual functions, improving their effectiveness, ensuring greater overall visual comfort and an increase in the benefits obtained with visual aids.

In some cases, Visual Training is necessary to improve the effectiveness of visual skills which could not be achieved with the sole use of glasses and/or contact lenses.



Exercises to improve eye coordination and binocular vision in usual activities with the brock string.



Exercises to improve peripheral vision and visual tracking with the Marsden sphere.



Exercises to improve low fusional and accommodative reserves and poor stereopsis with the Aperture rule.

Exercises to improve the effectiveness of visual skills


VISUAL TRAINING HELPS:

- To read better, faster, more effectively, without getting tired, reducing any tension in the muscles of the neck, shoulders and back in general, memorizing, understanding and therefore assimilating more easily in less time.
- To maintain a correct posture with greater balance and greater agility during any motor act, in particular during sporting activity, making it more precise, faster, effective and satisfying.
- To interpret visual information more effectively, creating a more adequate representation of space and the environment and therefore greater safety and reliability in driving vehicles.
- To support the complex organization of the extrinsic ocular muscles in the conduction of the most effective movements for reading during study and daily work in the office and in all activities with a high attention requirement, significantly reducing asthenopia.
- To synergistically integrate all the senses together (multisensory integration).

APPLICATIONS OF VISUAL TRAINING

The fields of application are varied and multiple:

EXAMPLES OF APPLICATIONS OF VISUAL TRAINING



Visual Training is also particularly effective in Global Postural Reprogramming that concerns the repercussions on the entire musculoskeletal system and "aims" to prevent or correct the postural anomaly.



Test to exercise accommodative and convergence flexibility

7.3 VISUAL TRAINING



Marsden spheres

TO IMPROVE THE ACCURACY AND SPEED OF TRACKING MOVEMENTS

New spheres with letters and numbers for visual training exercises useful for improving peripheral vision and visual tracking. Training with the Marsden sphere improves body coordination and perception of space. The exercises can be performed lying down, sitting, standing or on a platform in monocular or binocular vision. Equipped with a rope for hanging from the ceiling.

Ref. 08437 **Kr 471,-**

Soft sphere with letters on a white background for checking the ocular movement, fixation, accommodation and free space dynamics.

Dimensions: Ø 70 mm
Weight: 90 g
Rope length: 6 m
Material: foam sponge
1 pc.



Ref. 08438 **Kr 471,-**

Soft sphere with letters and numbers on a yellow background for eye motor activities, fixation, accommodation and free-space dynamics. The yellow background helps in perception of the letters with the contrast of the common white wall background.

Dimensions: Ø 70 mm
Weight: 90 g
Cord length: 6 m
Material: foam sponge
1 pc.



Ref. 08439 **Kr 551,-**


Soft sphere with letters on a Red/Green background for activities such as ocular movement, fixation, accommodation and free space dynamics hung from a string. To be used with the Red/Green goggle (Ref. 08291).

Dimensions: Ø 70 mm
Weight: 90 g
Rope length: 6 m
Material: foam sponge
1 pc.





Near Point Dot Cards



Ref. 08356 Kr 296,-

The purpose of this test is to measure the near point convergence. It can be used for home exercise. Very effective test as it ensures correct ocular alignment. Comprises of a narrow strip of double-sided paper with central fold. On one side there is a line of 10 points positioned at 3cm intervals along the central fold, on the other side there are two columns of letters 4 cm apart. Each pack contains 50 cards plus instructions.

Dimensions: 75x295 mm



MD MD Medical device compliant with EU Regulation no. 2017/745

Brock String



TO DEVELOP KINAESTHETIC AWARENESS OF CONVERGENCE/DIVERGENCE

Ref. 08416 Kr 601,-

The Brock String with 5 coloured beads is used for physiological diplopia training.

Used for jump vergences and other anomalies of binocular vision.

Dimensions:

Beads Ø 15 mm


String 4 m

6 pcs.

7.3 VISUAL TRAINING



Red/Green reading sheets



TO CHECK FOR THE PRESENCE OF OCULAR SUPPRESSION

Ref. 08275 **Kr 1582,-**

A Red/Green Reading sheets for anti-suppression training. The patient reads or looks at a picture wearing anaglyph glasses and occluding one eye at a time.

The sheets are to be held 2,5 cm above from the page or image. The sheets can be gently cleaned by rubbing them with a mild detergent or Isoproponal Alcohol or Centrostyle lens cleaning spray with a microfibre cloth. Do not use a paper towel as it will scratch the surface.

Dimensions: 127 x 178 mm

Bar width 15 mm

6 pcs.

NEW

Ref. 08275S **Kr 369,-**

Red and green individual reading sheet

1 pc.

Red and Green vynil sheets



TO CHECK FOR THE PRESENCE OF OCULAR SUPPRESSION

Ref. 08277 **Kr 1463,-**

Red and green sheets used for anti-suppression in visual training. Place sheets in front of tables, images and digital screens.

PRIMARY FILTERS

Green: Transmission Y%=3.6/Absorption 1.44

Red: Transmission Y%= 11.09/Absorption 0.063888

Dimensions: 216 x 280 mm

6 pairs

NEW

Ref. 08277G **Kr 248,-**

Individual green vynil sheet

1 pc.

NEW


Ref. 08277R **Kr 248,-**

Individual red vynil sheet

1 pc.

MD MD Medical device compliant
with EU Regulation no. 2017/745

Hendrickson lifesaver card



TO EXERCISE ACCOMMODATIVE FLEXIBILITY AND CONVERGENCE

Ref. 08278 **Kr 337,-**

White

6 pcs.

Ref. 08279 **Kr 478,-**

Transparent

6 pcs.

Two cards, one clear the other opaque.

Each containing two vertical line of four green lifesaver circles and four red lifesaver circles.

The user is required to demonstrate the maximum accommodative flexibility, convergence and accommodation-convergence interaction.

To develop speed, resistance and flexibility.

Prismatic glasses can be used with R/G and polarized filters.

Dimensions:107x158 mm

Red & green glasses for children and adults



GLASSES USED FOR ANAGLYPHIC TESTS

Ref. 08280 **Kr 948,-**

Children

1 pc.

Ref. 08291 **Kr 690,-**

Adult

1 pc.

NEW

The red/green anaglyphic glasses are ideal for use in several of the binocular vision tests as well as in vision control in red/green anti-suppression exercises.

The wrap-around function limits interfering side light. Ideal for most Lancaster Red-Green tests, Hess Screen Test colors and commonly used visual therapy, exercises requiring patients to wear red/green lenses.


Available in adult and children sizes.



MD MD Medical device compliant with EU Regulation no. 2017/745

7.3 VISUAL TRAINING

Transparent and opaque cards with eccentric red/green circles



NEW

Ref. 08272 Kr 344,-

Pairs of cards, transparent and opaque, with eccentric red/green circles to control central suppression and fixation disparity.

With the use of anaglyphic glasses, it is possible to carry out convergence and divergence exercises.

Dimensions: 50x76 mm

2 pcs. transparent cards

2 pcs. opaque cards.

Transparent and opaque cards with eccentric circles

NEW


Ref. 08274 Kr 326,-

Pairs of cards, transparent and opaque, with eccentric circles to control central suppression and fixation disparity.


Dimensions: 50x76 mm

2 pcs. transparent cards

2 pcs. opaque cards.



Red/blue dot cards for convergence exercises



NEW

Ref. 08276 Kr 134,-

Double-sided card with 3 red circles on one side and 3 blue circles on the other side.


Used to develop control and accuracy of vergences (convergence and divergence). Improve voluntary convergence and normalize the next point of convergence (PPC).

Useful in the presence of convergence insufficiency. It can be used to explain and improve the sense of physiological diplopia.

Dimensions: 63x140 mm

1 pc.

Binocular vision assessment test



VISION WITH LEFT EYE
RIGHT EYE COVERED



VISION WITH RIGHT EYE
LEFT EYE COVERED

FIXATION SUPPRESSION TARGET FOR ASSESSING BINOCULAR VISION AT DISTANCE

Ref. 08431 Kr 3766,-

- Quick test for assessing the presence of monocular suppression at distance, about 6 m, in free space or behind the phoropter.
 - Fixation target E in the centre of the test to test the ability to fuse images at a distance.
 - 4 arrows of two sizes to assess the presence of suppression in the right and left eye.
 - The test is performed in a well-lit area without glare.
- Dimensions: 145x127 mm

Complete with polarized spectacles (Ref. 08412.1) and instructions



VISION WITH LEFT EYE
RIGHT EYE COVERED



VISION WITH RIGHT EYE
LEFT EYE COVERED

FIXATION SUPPRESSION TARGET FOR ASSESSING BINOCULAR VISION AT NEAR

Ref. 08430 Kr 2959,-

- Quick test for assessing the presence of monocular suppression at close range, 41-46 cm, in free space or behind the phoropter.
 - E fixation target in the centre of the test to test the ability to fuse images at proximal distance
 - 4 arrows of two sizes to assess the presence of suppression in the right and left eye
 - The test is performed in a well-lit area without glare.
- Dimensions: 114x114 mm

Complete with polarized spectacles (Ref. 08412.1) and instructions



MD MD Medical device compliant with EU Regulation no. 2017/745

7.3 VISUAL TRAINING

Kit Aperture Rule



NEW

FUSIONAL TRAINER RULE

Ref. 08357 Kr 2394,-

In most visual training programs, the Fusional Trainer is a standard tool. This kit features progressive fixed demand and helps improve LOW fusional and accommodative reserves, poor stereopsis, fixation disparity suppression. No polarized or tinted lenses required.


The optometrist assists patients with eye strain and visually related learning difficulties.

Base-Out 30PD to 17.5PD Base-In.


The kit includes:

- Sport/Action Card Set
- Single and double opening cursors
- Red pointers
- Manual


EXAMPLES



Using the single opening for fusion reserve convergence exercises for exoforia



Using the double aperture for fusion reserve exercises of divergence for exoforia



Using the single opening for fusion reserve convergence exercises for exoforia

MD MD Medical device compliant with EU Regulation no. 2017/745

Bull's Eye



TO IMPROVE ACCOMMODATIVE FLEXIBILITY Ref. 08270 Kr 941-

Transparent sheet with concentric circles on the top and a matrix of letters and numbers at the bottom. The purpose of this target is to improve the eyes' focusing system. It helps build-up flexibility of accommodation quickly and accurately. Can help in sports performance skills to improve reaction time and hand-eye coordination and can also improve classroom reading and copying skills. The Bull's eye target must be held in front of the examined eye at a distance between 30 and 40 cm.

Dimensions: 70x146 mm
 10 pcs.

Hart's Tables



TO IMPROVE ACCOMMODATIVE FLEXIBILITY Ref. 08271 Kr 517,-

White
 Dimensions: 215 x 280 mm
 20 pcs.

Ref. 08273 Kr 633,-

Transparent
 Dimensions: 76 x 90 mm
 6 pcs.

Two tables, both with white background and black lettering, one is transparent. The purpose of this table is to help train the eye accommodative system from distance to near on demand. There are different ways to apply these charts depending on the result required (for example to develop the ability to perform rapid and accurate saccadic eye movements and to increase accommodative flexibility).

It can be used in conjunction with the Bull's Eye Test. The table set comprises of two cards of different size letterings.

MD MD Medical device compliant with EU Regulation no. 2017/745

7.3 VISUAL TRAINING

Wolff's double sphere target



STAINLESS STEEL

Ref. 08350 **Kr 259,-**

Two chromed fixation spheres of different diameters for quick evaluation of ocular motility, fixation stability, saccadic movements and pursuits.

Used in the NPC test (near point of convergence for the evaluation of the break and recovery of convergence).

They can use by the patient themselves as a fixation target.

Spheres: Ø 6 and 13 mm - stick: 36 cm
2 pcs.



NEW

RED/GREEN

Ref. 08351 **Kr 403,-**

One set of RED and GREEN fixation spheres of different diameters. Used for quick assessment of eye movement, fixation stability, saccadic movements and pursuit movements (pursuits with anti-suppressive control in the NSUCO test procedure).

Spheres: Ø 6 and 13 mm - stick: 36 cm
2 pcs.



Lang Fixation Cubes - Lang Fixation Stick

LANG FIXATION CUBES

Fixation cube white and red with pictures designed to attract the attention of children in order to observe the eyes motility and the fixation movements in all directions of gaze. They are being used for assessment of eye fixation, motility, accommodation, convergence, for cover test and to prepare the patient for testing with the LANG-STEREOTEST. Usage at 30 cm reading distance.

Ref. 08352R

Kr 483,-

Red

Length 155 mm - 1 pc.

Ref. 08352W

Kr 483,-

White

Length 155 mm - 1 pc.



LANG FIXATION STICK

Ref. 08355 **Kr 383,-**

FRONT: Fixation Test with the same pictures as on the Lang Fixation Cubes

BACK: Examination of accommodation and convergence with numbers.

Recognizing the pictures prepares the child for the Lang I and II tests.

Visus: Distance of 30 cm 0.2-0.5.

Length 146 mm - 1 pc.



MD Medical device compliant with EU Regulation no. 2017/745

7.3 VISUAL TRAINING

Matt black/Frosted Flip-Up occluder glasses

Children's flip-up occluder glasses are made of TPE (thermoplastic rubber) and latex-free. Each occluder flips up and down for easy occlusion with large circular lenses for maximum coverage. The new gender neutral colours make them suitable for boys and girls. The rubber temple tips help prevent the glasses from slipping. The flexible and durable plastic makes these glasses almost indestructible. Available in matt or frosted occluder lenses.



Ref. 08282 **Kr 832,-**
Frosted
1 pc.



Ref. 08283 **Kr 832,-**
Matt black
1 pc.

Glasses with R/G filter - Tiger shape

NEW

Ref. 08287 Kr 1076,-

Nice children's glasses in soft material with red/green filters for carrying out anaglyphic tests.

From 4 years old.

1 pc.



Glasses with polarized lenses - Parrot shape



NEW

Ref. 08288 Kr 1076,-

Nice children's glasses in soft material with polarized lenses to be able to carry out various stereo tests.

From 4 years old.

1 pc.

MD MD Medical device compliant
with EU Regulation no. 2017/745

7.3 VISUAL TRAINING



Home Trainer Kit




Ref. 08285 **Kr 2052,-**
Tube designed to hold a selection of 16 home training products.

Kit includes:

- 1 Bull's eye target;
- 3 accommodative charts with letters, numbers and arrows (9, 18, 72 point chart), made of cardboard;
- 3 antisuppressive chart with letters, numbers, and arrows (9, 36, 72 point chart), made of cardboard;
- 1 red-green goggle;
- 1 couple of fusion cards;
- 1 peripheral awareness chart with letters, made of cardboard;
- 1 red-green bar reader;
- 1 soft training ball with letters, adhesive hook to hang the ball enclosed;
- 1 trainer string;
- 1 couple of triangular pencil trainers;
- 1 translucent polyfusion;
- 1 translucent occluder.

Everything is contained in a sturdy cardboard tube.
Tube dimensions: Ø 11 cm., length 43 cm., weight 600 gr.

How to use the Home Trainer Kit



HART'S TABLETS - BULL'S EYE


OCULOMOTOR PROCEDURES

Fixation - Rotation - Pursuit

They improve the motor skills of the visual system by organising space and pursuing a target.

ACCOMMODATIVE PROCEDURES

To **increase the fusional power** and **flexibility of accommodation** and build control and awareness of the accommodative system.



ANAGLYPHIC TESTS

3 RED-GREEN POLYHEDRAL BOARDS

1 ANAGLYPHIC GOGGLE

Useful for **developing the ability to move the focus point quickly**, from near to far, while **maintaining good eye coordination without diplopia and suppression**.

The tests are carried out using the anaglyphic goggle.



ANAGLYPHIC TESTS

2 FUSION EXERCISE CARDS

1 BOARD WITH LIFESAVER RINGS

1 BOARD WITH RED-GREEN BARS

Useful for **developing the ability to use the eyes in a coordinated way** and **achieve flexible interaction** between **the focusing system** and **the eye coordination system**.

The tests are carried out with the use of anaglyphic goggles.



BINOCULAR EXERCISE TEST

1 MARSDEN TRAINING BALL

To **develop** awareness of spatial assessment and **focusing skills** used to improve accuracy and speed of tracking and turning movements.

2 PENCILS FOR TRAINING EXERCISES

For **improving coordinated eye function** over long periods.

1 BROCK STRING

For **improved eye coordination** and increased ability to **move binocular vision** from one point to another in space quickly and accurately.

1 OCCLUDER - EYE PATCH WITH ELASTIC CORD

For use in **monocular assessments**.

7.3 VISUAL TRAINING



Lectern for a correct posture



Ref. 08286 **Kr 1440,-**

The lectern has the following features:

- Supports A4 and A3 sheets and books
- Smooth and rigid surface for writing and drawing
- Designed not to hinder the movement but to support the forearm in writing for both for right-and left-handed people,
- Its low slope support minimises sheets, book and/or pencils slippage
- A handle for easy carrying
- The upper part shows the correct positions to adopt and the incorrect ones to avoid; in addition, these drawings represent peripheral visual stimulation
- The supports feet are easily assembled by removing the protective film from the double-sided adhesive that is already applied.

The kit includes:

- ergonomic grip, pencil, functional instructions for the user, assembly instructions.

Top dimensions: 485 x 370 mm

Inclination: 15°


Total weight: 500 g





Lectern for a correct posture

The inclined plane is suitable for use at school and home to ensure support for the correct angle during visual-cognitive exercise and to impose correct posture placement.



THE RULE OF THE NUMBER 3

The pen must be kept between thumb, index and middle index that acts as a support. The fingers form an equilateral triangle.

ERGONOMIC GRIP:

Used for fixation exercises. The ergonomic triangular silicone sleeve fits over pens or pencils for comfortable grip by the hand to support correct posture when writing. The triangular shape of the handle helps to maintain the correct grip of the pencil during drawing and writing activities. Has black alphabetic characters printed on the body pencil and the ergonomic grip.

Ref. 08286R/B Kr 59,-

Pack:

2 pcs. Ergonomic grips - BLUE

2 pcs. Ergonomic grips - RED



7.4 THE BEGINNING OF VISUAL ANALYSIS

Functional tests

The evaluation of the functional state of the subject consists of a series of tests that investigate the different abilities of the visual system in recognizing objects around it, such as visual acuity, contrast sensitivity and spatial localization. In order to maintain clear and comfortable binocular vision at all viewing distances, each individual requires a number of well-functioning skills:

- must be able to align the two eyes and maintain alignment for a prolonged period of time;
- must have accurate and efficient accommodation based on the request and must be able to sustain it comfortably over time;
- accommodation and convergence must interact appropriately.

All this cannot ignore a state of good health of the visual organs.

Given the different visual abilities, it is necessary to adopt various types of tests, some of which provide a predominantly quantitative evaluation of visual performance, while others emphasize the qualitative aspect.

To obtain maximum information from visual skills tests, it is important to take note not only of the quantitative result of the performance produced but also of the quality of the response.

The evaluation of the visual system can be divided into three main moments:

- 1 **Preliminary tests:** they aim to provide mainly qualitative information about the degree of development of binocular integration and the degree of functioning of the sensorial, integrative motor systems.



Find the prismatic value that compensates for ocular deviation using the orange occluder (Ref. 08348).



Test for ocular dominance



Convergence and divergence exercises




Using the trial frame

Functional tests

2 **Phoropter tests:** allow you to evaluate the subjective refraction of the eye. Through this instrument, it is possible to perform numerous visual tests thanks to the vast combination of auxiliary lenses that facilitate the realization of all possible refractive options that cannot be achieved through other analysis methods.



3 **Analysis in free space:** allows you to better observe the attitudes that the subject adopts during the analysis by providing qualitative information, reducing the sense of "claustrophobia" given by the phoropter, allowing a more natural posture of the body and eyes, as well as allowing a considerable increase in the peripheral field, more similar to the subject's real visual needs.



Fixations



Saccades




Pursuits


7.4 DIAGNOSTICS

Berens Prism Bar - Vertical and horizontal bar

VERTICAL BAR



HORIZONTAL BAR



The prism dioptres are used to evaluate strabismus. With the use of cover testing it is possible to measure, in prism dioptres, phorias and tropias.

Ref. 08300 Kr 3556,-

In a rigid protective case, complete with:

- 1 vertical bar with the following prism dioptres:
1^Δ, 2^Δ, 3^Δ, 4^Δ, 5^Δ, 6^Δ, 8^Δ, 10^Δ, 12^Δ, 14^Δ, 16^Δ, 18^Δ, 20^Δ, 25^Δ.
- 1 horizontal bar with the following prism dioptres:
1^Δ, 2^Δ, 4^Δ, 6^Δ, 8^Δ, 10^Δ, 12^Δ, 14^Δ, 16^Δ, 18^Δ, 20^Δ, 25^Δ, 30^Δ, 35^Δ, 40^Δ.



MD MD Medical device compliant
with EU Regulation no. 2017/745

Retinoscopy Racks



For dynamic retinoscopy use in conjunction with the Wolff's Target. The presence of positive and negative lenses allows to neutralize both main meridians. If used in conjunction with a phoropter the field of view can be reduced by 70%.

Ref. 08304 Kr 6429,-

In a rigid protective case, complete with:

- 1 red rack from -1.0D to -10.0D with additional sliding lens -0.5D and -10.0D
- 1 green rack from +1.0D to +10.0D with additional sliding lens +0.5D and +10.0D



MD MD Medical device compliant with EU Regulation no. 2017/745

7.4 DIAGNOSTICS

Plates used for determining eye fixation dominance

MATT WHITE WITH MIRROR



TRANSPARENT

NEW

Ref. 08346 Kr 175,-

Two plexiglass plates, transparent and matt, with a hole in the center (approx. 2 cm).


They are used to assess which eye is the dominant eye in monocular fixation of a sight, from a distance and at the proximal distance using the mirror strip.

It helps identify the eye on which to preferably place the filters and dissociating prisms during the execution of the preliminary tests.

NOTE: The transparency of the plexiglass allows the maintenance of the fixation of the aim at a distance while raising the arms in order to center the aim through the hole.


On the other hand, a matt plate, does not allow observation of the aim while raising the arms, thus increasing the probability of error in centering the examining eye.

Dimensions:
253x153 mm
2 pcs.



Instructions available in Italian, English, French, German, Spanish.

Ishihara Colour Test for colour deficiency 10-14-24-38 plates



ISHIHARA COLOUR TESTS

38 PLATES

Ref. 08400 **Kr 7384,-**

Used principally by ophthalmologist.

24 PLATES

Ref. 08401 **Kr 5472,-**

Used principally for occupational screening.

14 PLATES



Ref. 08399 **Kr 3009,-**

Used mainly for driving licence renewals.

10 PLATES

Ref. 08402 **Kr 3625,-**

Used principally for children and illiterates.



The Ishihara colour tests are:

- Internationally recognised
- Pseudoisochromatic plates
- Simple to use
- Speedy to use
- Efficient in detecting red/green colour vision deficiency
- Includes instruction booklet

To prevent colour fading and to keep hues sharp never expose plates to bright sunlight or leave pages open for any length of time. Test should be conducted with light levels at 500-600 lux.



MD MD Medical device compliant with EU Regulation no. 2017/745

7.4 DIAGNOSTICS

Binocularity degrees

The development of Binocularity occurs gradually, and it is possible to divide it into different steps. Worth distinguished three phases and named them "degrees of development of binocular vision":

1° degree - **Simultaneous perception or Biocularity**: the simultaneous perception of the two eyes without the possibility of binocular fusion. Through this phase it is possible to investigate the presence or absence of suspension or suppression and to detect heterophoria. Test example: Test with prism Von Graefe.

2° degree - **Flat fusion or Binocularity**: two images strike the two retinas in corresponding points, after

having reached the occipital visual cortex, and are associated in a single image that represents the fusion of the two primitive images. Test examples: Worth's lights or Red Filter Test.

3° degree - **Stereopsis**: a visual-perceptive process which, starting from two slightly different monocular retinal images, leads to the sensation of depth. It occurs only at the level of the striated cortex where the re-elaboration of the information of the two eyes occurs. Examples of stereotests; **Lang, Titmus, Butterfly**.

1° degree - Simultaneous perception or Biocularity



PRISM 6 Δ

Ref. 08307 **Kr 791,-**

Prism of 6 prismatic diopters (Δ) useful to detect the first degree of biocularity.

The prism must be placed with a vertical base in front of one of the eyes while the person is fixing a point-like target projected from a distance. A prism

of this power cannot be compensated for by vertical fusion reserves, therefore, under normal conditions, the person should report vertical diplopia. Otherwise it could be a case of suppression. Before proceeding with further investigations, move the prism to the other eye or flip the base

1 pc.

MD MD Medical device compliant
with EU Regulation no. 2017/745

2° degree - Flat fusion or Binocularity

RED FILTER - PENLIGHT - OCCLUDER

The red occluder determine a level of dissociation between the two eyes. The orange filter, less dissociating than the red one, helps the practitioner to more adequately quantify the prescriptive value for restoring single binocular vision. The purpose of the test is to evaluate the fusional stability (hence the quality) of

binocular vision at different distances and in different gaze positions. The red filter test can highlight the presence or absence of suspension or suppression of the monolateral and alternating type. It is recommended to use the penlight and the black occluder to carry out the test.



Ref. 08311
Kr 98,-
Red
occluder 1 pc.



Ref. 08349
Kr 98,-
Shiny black
occluder 1 pc.

Ref. 08317
Kr 127,-
White
penlight 135
mm - 1 pc. 2
batteries AAA
1,5V included.



WORTH 4 DOT FLASHLIGHT WITH RED-GREEN GLASSES

Ref. 08305 **Kr 1623,-**
Flashlight with 4 coloured dots (one red, two green and one white) with a diameter of 6 mm each and anaglyph glasses with red-green filters. Test performed to assess the presence of flat binocular fusion and its stability in different gaze positions from near and to detect the presence or absence of

monolateral and/or alternating suspension or suppression. The test can be performed while wearing prescription lenses or with trial lenses, in low lighting conditions. Worth's test should also be performed in cases of reduced visual acuity that does not improve with the pinhole test.



Flashlight battery type D
not included.

MD MD Medical device compliant
with EU Regulation no. 2017/745

7.4 DIAGNOSTICS

Hand held occluders



Ref. 08308 **Kr 574,-**

Set of 5 occluders with plexiglass display:

- 1 red occluder,
- 1 frosted occluder,
- 1 black paddle occluder,
- 1 multi pinhole occluder,
- 1 black shiny occluder.

Dimensions: Display: 245x120 mm
Occluders: Ø 60 mm x 240 mm



Ref. 08349
Kr 98,-
Shiny black occluder



Ref. 08314
Kr 109,-
Multi pinhole
occluder Holes
Ø 1.4 mm



Ref. 08313
Kr 98,-
Black occluder



Ref. 08312
Kr 116,-
Frosted occluder



Ref. 08311
Kr 98,-
Red occluder

Hand held occluders



NEW

Ref. 08348 **Kr 95,-**
Orange occluder

In situations of diplopia, the search for the prismatic value that compensates for the ocular deviation in the considered gaze position should be sought by minimizing the dissociation between the two eyes. The orange filter was created for this purpose because it is less dissociating than the classic red filter used to identify the maximum deviation as it has a high chromatic dissociation. Using the orange filter can help the clinician to more adequately quantify the prescribable value for the restoration of single binocular vision since, as is well known, the more the degree of dissociation between the two eyes increases, the greater the angle of deviation found will be.



NEW

Ref. 08347 **Kr 95,-**
Matt white occluder

Used so as not to alter the perception of ambient light intensity and not to create pupil dilation. Occlusion is however total and the test appears more ecological.

Red Maddox Rod



Ref. 08316 **Kr 547,-**

Red Maddox Rod with grooves at 45° from the centre of the lens.

Used in combination with prisms to measure horizontal and vertical phorias from near and from far.

Length: 220 mm

Lens diameter: 50 mm



MD Medical device compliant
with EU Regulation no. 2017/745

7.4 DIAGNOSTICS


3° degree - Stereopsis

FLY STEREOPSIS TEST

Ref. 08412 Kr 4788,-

- Rapid test for amblyopia and strabismus
- Traditional Fly picture testing Gross Stereopsis
- 10 levels of 4 circle disparity for critical testing
- Graded circle test from 400 seconds down to 20 seconds (with no monocular clues)
- Fine stereopsis (400 to 100 seconds of arc) measuring 3 Gross levels of disparity using the internationally recognised Symbols - Square, Circle, Plus and Triangle
- Ideal for both children and adults including non-reading and non-verbal
- New easy-to-hold booklet
- Answer key on back cover
- Instructions included.

Includes 1 adult and 1 child's polarized goggles




STEREO BUTTERFLY TEST

Ref. 08413 Kr 4900,-

- Rapid test for amblyopia and strabismus
- Traditional Butterfly picture testing Gross Stereopsis
- 10 levels of 4 circle disparity for critical testing
- Graded circle test from 400 seconds down to 20 seconds (with no monocular clues)
- Fine stereopsis (400 to 100 seconds of arc) measuring 3 Gross levels of disparity using the internationally recognised Symbols - Square, Circle, Plus and Triangle
- Ideal for both children and adults including non-reading and non-verbal
- New easy-to-hold booklet
- Answer key on back cover
- Instructions included.

Includes 1 adult and 1 child's polarized goggles



MD MD Medical device compliant with EU Regulation no. 2017/745


3° degree - Stereopsis

RANDOM DOT STEREOPSIS TEST

Ref. 08414 Kr 4900-

- Rapid test for amblyopia and strabismus
- Ideal for children and adults who are non-readers and non-verbal
- Expanded Random Dot Symbols Test - Square, Circle, Plus and Triangle (400, 200, 100 seconds of arc).
- Graded circle test to 12.5 seconds of arc with no monocular clues
- New improved booklet
- Answer key on back cover
- Instructions included.

Includes 1 adult and 1 child's polarized goggles



Ref. 08412.2

Kr 766,-

Children

1 pc.



Ref. 08412.1

Kr 533,-

Adult

1 pc.



MD MD Medical device compliant with EU Regulation no. 2017/745

7.4 DIAGNOSTICS

Lang Stereo Test® I-R and II-R - For children - Without the use of polarized glasses

The LANG STEREO TEST® is a rapid test for the screening of stereoscopic vision disorders in children. This test provides early diagnosis without the use of polarized glasses.


LANG STEREO TEST® I-R

Ref. 08393 Kr 4700,

Three random point stereo objects with different disparities CAT 1200", STAR 600" and CAR 400". Test distance 30-40 cm Dimensions: 14.5x9.5 cm with rounded corners.

New features of the LANG STEREO TEST® I-R:

- INCREASED SHARPNESS AND UNIFORMITY OF THE PATTERN a 400" disparity, making the gradation between objects more harmonious.
- REDUCED SINGLE EYE RECOGNITION. In monocular vision these stereograms do not reveal any contours, while in binocular vision the areas that produce a horizontal disparity are seen in relief.
- ROUNDED CORNERS
- BRIEF INSTRUCTIONS IN ENGLISH ON THE BACK OF THE TEST CARD, WITH SCHEMATIC DIAGRAM FOR USE.



LANG STEREO TEST® II-R

Ref. 08394 Kr 4700,-

Four 3D transverse disparity test figures ELEPHANT 600", JEEP 400", MOON 200" and STAR 200" (always visible). Test distance 30-40 cm.

Size: 14.5 cm x 9.5 cm with rounded corners.

New features of the LANG STEREO TEST® II-R

- IMPROVED OPTICAL AND PRINT QUALITY
- THREE FIGURES VISIBLE ONLY BINOCULARLY (moon, jeep and elephant) AND ONE SINGLE VISIBLE FIGURE (star). The disparities are chosen to be finer
- ROUNDED CORNERS
- BRIEF INSTRUCTIONS IN ENGLISH ON THE BACK OF THE TEST SHEET, WITH SCHEMATIC DIAGRAM FOR CORRECT USE (specular arrangement of stereophonic objects and typical pattern of eye movements in stereophonic subjects).



MD MD Medical device compliant with EU Regulation no. 2017/745

Teddy Bear Test - 3° degree - Stereopsis



TEDDY BEAR TEST

Ref. 08417 **Kr 4900,-**

Accurate and fast test for the determination of amblyopia and strabismus. Also used for convergence and divergence exercises.

- Teddy Bear Vectograph has a range of 2,600 seconds of arc stereopsis depth.
- Teddy Bear & Blocks orientation at +30 to +1500 seconds of arc and -25 to -1100 seconds of arc.
- Used to strengthen the binocularity system and to provide Base-In & Base-Out training
- 16-diopter range for base-in (divergence) training and 24-diopter range for base-out (convergence) training for a total accommodation range of 40 diopters.
- Set includes: Protective binder with pen, Doctor and Patient Manuals, Therapy Record Form & Standard Polarized Goggle.

Complete with polarized spectacles for adult
(Ref. 08412.1)

1 pc.

7.4 DIAGNOSTICS


Amsler Test

centrostyle®
around PEOPLE

ISTRUZIONI DI UTILIZZO


1. Se richiesto, fare questo test con gli occhiali da lettura indossati tenendo il foglio a 30 cm.

2. Coprite un occhio con una mano e fissate il punto bianco al centro del test per qualche secondo.



3. Eseguite il test sull'altro occhio.

4. Se notate che qualche linea o quadrato appaiono ondulati, distorti, irregolari o se durante il test dovete percepire delle zone scure, consultate appena possibile uno specialista.



centrostyle®
around PEOPLE

INSTRUCTIONS

1. If necessary, wear reading glasses and hold the card at 11 in (28 cm).



2. Cover one eye with your hand and look at the white square in the center. If the lines or squares appear wavy, distorted, irregular or if you notice dark areas during the testing of one eye, consult a licensed eye professional.

ANWEISUNG

1. Wenn notwendig, tragen Sie Lesegläser und halten Sie das Blatt 30 cm von dem Auge.



2. Bedecken Sie ein Auge mit einer Hand und schauen Sie auf den Weißpunkt in der Mitte für einige Sekunden.

3. Wenn Sie bemerken, dass einige Zeilen wellig, deformiert, unregelmäßig oder dunkle Zonen erscheinen, sollten Sie man dunklen Stellen bemerken, bitte wenden Sie sich an einen Spezialisten.

INSTRUCCIONES DE USO

1. Si necesario, hacer la prueba con gafas de lectura y sosteniendo la hoja a 30 cm.



2. Cubrir un ojo con una mano y fijar el punto blanco en el centro de la prueba durante unos segundos.

3. Hacer la prueba de la misma forma en el otro ojo.

4. Si durante la prueba alguna línea o el cuadrado parecen ondulados, distorsionados, irregulares o si percibe áreas oscuras, consulte a un especialista.

MODE D'EMPLOI

1. Si nécessaire, effectuez ce test avec lunettes de lecture portées tenant la feuille à 30 cm.



2. Couvrez un œil d'une main et fixez le point blanc au centre du test pendant quelques secondes.

3. Faites le test sur l'autre œil.

4. Si vous remarquez que quelques lignes ou carré semblent ondulés, déformés, irréguliers ou si pendant le test vous devriez percevoir des zones sombre, consultez un spécialiste.

Instructions available in Italian, English, French, German, Spanish.

TO QUALITATIVELY TEST THE CENTRAL VISUAL FIELD

Ref. 08409 Kr 95,-


Heavy card consisting of a white grid on a black background to check changes in the macula.

It is a qualitative, fast-performing technique to test the central visual field (10°) for the presence of age-related macular degeneration (ARMD) by identifying a relative or absolute scotoma.

The subject is asked to describe what he sees, if the lines distort or disappear, while the fixation remains stable on the central point.

The test must be performed without wearing bifocal or multifocal lenses, occluding one eye at a time and with natural lighting.

Dimensions: 130 x 175 mm
2 pcs.



Binocular confirmation vision test - Cross cylinders



Binocular confirmation test

BINOCULAR CONFIRMATION TEST

Ref. 08325 Kr 1504,-

Used to verify that the dioptré previously tested on both eyes is adequate.

With two sliding lenses for optical centring, in protective case.

- Powers: $\pm 0.25D$
- Powers: $\pm 0.50D$



Cross cylinders

CROSS CYLINDERS

Ref. 08330 Kr 574,-

Used to refine axis and power of astigmatism cylinders. With knurled handles in protective case.

- Powers: $\pm 0.25D$
- Powers: $\pm 0.50D$

7.5 PLASTIC FLIPPERS

Binocular visus confirmation test - Optical crystal lenses

 **Ref. 08331** Kr 109,-
Flipper to place lenses for test purposes

 **Ref. 08332** Kr 378,-
Confirmation Test $\pm 0.25D$

 **Ref. 08333** Kr 378,-
Confirmation Test $\pm 0.50D$

 **Ref. 08334** Kr 378,-
Confirmation Test $\pm 0.75D$

 **Ref. 08335** Kr 378,-
Confirmation Test $\pm 1.00D$

 **Ref. 08336** Kr 378,-
Confirmation Test $\pm 1.50D$

 **Ref. 08337** Kr 378,-
Confirmation Test $\pm 2.00D$

 **Ref. 08341** Kr 378,-
Confirmation Test $\pm 2.50D$

 **Ref. 08342** Kr 378,-
Confirmation Test $\pm 3.00D$

 **Ref. 08338** Kr 492,-
Flipper with polarized lenses and Red-Green filters

Ref. 08339 Kr 2724,-
Set of 6 flippers with display
• 5 flippers with sphere lenses in optical glass.
Powers: $\pm 0.25D$, $\pm 0.50D$, $\pm 1.00D$, $\pm 1.50D$, $\pm 2.00D$
• 1 flipper with polarized lenses and Red-Green filters.

Ref. 08329 Kr 3513,-
Set of 9 flippers with display
• 8 flippers with sphere lenses in optical glass.
Powers: $\pm 0.25D$, $\pm 0.50D$, $\pm 0.75D$, $\pm 1.00D$, $\pm 1.50D$,
 $\pm 2.00D$, $\pm 2.50D$, $\pm 3.00D$
• 1 flipper with polarized lenses and Red-Green filters.



PRISMATIC FLIPPERS

7.5

Plastic flippers with optical crystal lenses

Mainly used for exercises on the flexibility of convergence/divergence. As the prismatic power increases so does the difficulty of the exercise

Used in the test of binocular balance provided for by the 21-point procedure.

Used especially in cases of anisometropia to balance the optical compensation detected with the subjective examination monocular distance view.



NEW

Ref. 08322 Kr 752,-

Prismatic flipper 4BI $\blacktriangleright\blacktriangleleft$ / 6BO $\blacktriangleleft\blacktriangleright$
 4 BI means (2 BI+2 BI on one side).
 6 BO means (3 BO+3 BO on the other side).



NEW

Ref. 08327 Kr 752,-

Prismatic flipper 3BU \blacktriangledown +3BD \blacktriangle / 3BD \blacktriangle +3BU \blacktriangledown



NEW

Ref. 08323 Kr 752,-

Prismatic flipper 8BI $\blacktriangleright\blacktriangleleft$ / 12BO $\blacktriangleleft\blacktriangleright$
 8 BI means (4 BI+4 BI on one side).
 12 BO means (6 BO+6 BO on the other side).



08322



08323



NEW

Ref. 08324 Kr 752,-

Prismatic flipper 10BI $\blacktriangleright\blacktriangleleft$ / 10BO $\blacktriangleleft\blacktriangleright$
 10 BI means (5 BI+5 BI on one side).
 10 BO means (5 BO+5 BO on the other side).



08324



08327

MD MD Medical device compliant with EU Regulation no. 2017/745

7.6 CHILDREN'S TRIAL FRAME

0-3 years - For optometrical refractions - Tabo System

□ 34) 15 — 110

- PD: 50 mm
- Weight: 16 gr.


Adjustable Silicone Tube



Three lens holders suitable for standard ø 38 mm trial set lenses.

Wrap-around temples

Frame without any metal parts



Ref. 15731 Kr 307,-
Light Blue



Ref. 15730 Kr 383,-
Pink

• **COMFORTABLE AND LIGHTWEIGHT**

The design of the frame and its hypoallergenic material have been manufactured to offer the best possible visual comfort during eye examination. The wrap-around temples give a good fit to the child's face.

• **OPTICAL QUALITY**

The bridge is centred and splayed, thereby helping to maintain the optical centres. The frame has three external lens holders suitable for all ø 38 mm trial set lenses.

Determination of the astigmatism axis is according to the **TABO system**.

MD MD Medical device compliant with EU Regulation no. 2017/745



CHILDREN'S TRIAL FRAME

7.6

Ultra-light - For optometrical refractions - Tabo System



CHILDREN'S TRIAL FRAME KIT

Ref. 08215 Kr 1860,-

Display box with 4 children's trial frames. PD range 50-52-54-56 mm

- Lightweight, non-allergic structure, suitable for all standard Ø 38 mm trial lenses.
- Standard axis rotation.
- Comfort bridge.
- Flex temples, length adjustable from 90 to 130 mm.
- 4 lens holders each eye: 3 external rotating and 1 fixed internal.
- Manual lens axis rotation, scale in 5° increments.



Weight: 30 g

MD MD Medical device compliant with EU Regulation no. 2017/745

7.6 TRIAL FRAME

Ultralight - For optometrical refractions- Tabo System

- 1 - Pupillary distance adjustment
- 2 - Left temple pantoscopic adjustment
- 3 - Left lens axis adjustment
- 4 - Bridge height adjustment
- 5 - Right lens axis adjustment
- 6 - Right temple pantoscopic adjustment



Ref. 08220 Kr 1860,-

Non-allergic, light weight and strong, for optometrical refractions. Suitable for all standard \varnothing 38 mm trial lenses. Standard axis rotation.

- Determination of the astigmatism axis is according to the TABO system.
- Temples adjustable in length and inclination.
- PD adjustment from 54 to 70 mm with selector in 2 mm increments.



- Adjustable pads.
- Bridge height adjustable (+4 mm).
- Temple length adjustable from 95 to 135 mm.
- Temple pantoscopic $\pm 10^\circ$ from the horizontal by means of a screw.
- 4 lens holders each eye: 3 external rotating and 1 fixed internal.
- Manual lens axis rotation, scale in 5° increments

Weight: 49 g

MD MD Medical device compliant with EU Regulation no. 2017/745

For optometrical refractions - Tabo System

- 1 - PD adjustment, right and left eye
- 2 - Bridge height adjustment
- 3 - Corneal vertex distance adjustment
- 4 - PD measurement scale
- 5 - Scale for reading the axis position and base position
- 6 - Scale for measuring the corneal vertex distance (VD)
- 7 - Temple inclination adjustment
- 8 - Adjustable temple ends
- 9 - Temple length slide adjustment
- 10 - Knob for setting the axis position and base position
- 11 - Lens holder
- 12 - Nose pad on rack



Ref. 08228 Kr 4076,-

Lightweight structure, suitable for all standard Ø 38 mm trial lenses.

- Standard axis rotation
- Back vertex distance measurement gauge
- Independent monocular adjustments for optical centres
- Vertical and projection bridge adjustments
- Temple length and inclination adjustable
- PD adjustable range 50 to 80 mm
- Right and left PD adjustable independently from 25 to 40 mm
- Bridge height adjustable by up to 23 mm
- Temple adjustable in length from 100 to 135 mm and in inclination from +6° to -15°
- 5 lens holders each eye: 3 external rotating by 360° and 2 fixed internal
- Axis rotation is continuously adjustable, in increments of 5°.

Weight: 55 g


Ref. 08228.1 Standard nose pad - 1 pc.

Ref. 08228.2 Small nose pad - 1 pc.

Ref. 08228.3 Complete right temple - 1 pc.

Ref. 08228.4 Complete left temple - 1 pc.

Ref. 08228.5 Metal bridge - 1 pc.





MD Medical device compliant with EU Regulation no. 2017/745

7.6 TRIAL FRAME

For optometrical refractions - Tabo System - Made in Japan

- 1 - PD adjustment, right eye
- 2 - Bridge height adjustment
- 3 - PD adjustment, left eye
- 4 - Left temple inclination adjustment
- 5 - Left lens axis adjustment
- 6 - Vertical height adjustment
- 7 - Right lens axis adjustment
- 8 - Right temple inclination adjustment
- 9 - Back Vertex distance gauge (both sides)
- 10 - Temple length slide adjustment



LIGHTWEIGHT TRIAL FRAME

Ref. 08225 **Kr 4446,-**

Lightweight, non-allergic structure, suitable for all standard Ø 38 mm trial lenses.

- Standard axis rotation.
- Back vertex distance measurement gauge.
- Independent monocular adjustments for optical centres.
- Vertical and projection bridge adjustments.
- Temple length and inclination adjustable.
- PD adjustable range 48 to 80 mm.
- Right and left PD adjustable independently from 25 to 40 mm.
- Bridge height adjustable by up to 15 mm.
- Temples adjustable in length from 98 to 135 mm.
- Temple inclination adjustable from +6° to -15°.
- 5 lens holders each eye: 3 external rotating by 360° and 2 fixed internal.
- Axis rotation is continuously adjustable, in increments of 5°.

Weight: 58 g

Ref. 08225.1 **Kr 875,-**

Complete bridge - 1 pc.



MD MD Medical device compliant with EU Regulation no. 2017/745

TRIAL FRAME

7.6

For optometrical refractions - Tabo System

- 1 - PD adjustment, right eye
- 2 - Bridge height adjustment
- 3 - PD adjustment, left eye
- 4 - Left temple inclination adjustment
- 5 - Left lens axis adjustment
- 6 - Locking screw - left axis adjustment
- 7 - Bridge projection adjustment
- 8 - Locking screw - right axis adjustment
- 9 - Right lens axis rotation
- 10 - Right temple inclination adjustment



Ref. 08224 **Kr 1887,-**

Lightweight, non-allergic structure, suitable for all standard Ø 38 mm trial lenses. Standard axis rotation. Independent monocular adjustments for optical centres. Vertical and projection bridge adjustments for corrected vertex distances. Temple length and inclination adjustable.

Weight: 72 g


- PD adjustable from 48 to 80 mm
- Right and left PD adjustable independently from 24 to 40 mm.
- Bridge height adjustable by up to 13 mm.
- Bridge projection distance adjustable by up to 16 mm.
- Temples adjustable in length from 98 to 135 mm
- Temple inclination adjustment from +6° to -20°.
- 4 lens holders: 3 external and 1 internal.
- External axis rotation is continuously adjustable, scale at 5° increments.

MD MD Medical device compliant with EU Regulation no. 2017/745

7.6 TRIAL FRAME OCULUS UB4

For optometrical refractions - Tabo System

- 1 - PD adjustment, right and left eye
- 2 - Bridge height adjustment
- 3 - Corneal vertex distance adjustment
- 4 - PD measurement scale
- 5 - Scale for reading the axis position and base position
- 6 - Scale for measuring the corneal vertex distance (VD)
- 7 - Temple inclination adjustment
- 8 - Adjustable temple ends
- 9 - Temple length slide adjustment
- 10 - Knob for setting the axis position and base position
- 11 - Lens holder
- 12 - Polarizing filter
- 13 - Nose pad on rack
- 14 - Knob for turning the polarizing filter for 45° and 135°



Ref. 08234 Kr 10 602-

Lightweight, non-allergic structure, suitable for all standard Ø 38 mm trial lenses.

- Standard axis rotation
- Back vertex distance measurement gauge
- Independent monocular adjustments for optical centres
- Vertical and projection bridge adjustments
- Temple length and inclination adjustable
- PD adjustable range 48 to 80 mm
- Right and left PD adjustable independently from 24 to 40 mm
- Bridge height adjustable by up to 15 mm
- Temple adjustable in length from 98 to 135 mm and in inclination from +6° to -15°
- 5 lens holders each eye: 3 external rotating by 360° and 2 fixed internal
- Axis rotation is continuously adjustable, in increments of 5°.


Weight: 63 g (without polarizing filters)

Ref. 08234.1B Blue nose pad - 1 pc.

Ref. 08234.1G Grey nose pad - 1 pc.

Ref. 08235 Linear polarizing filter (1 pair)

Ref. 08236 Circular polarizing filter (1 pair)




MD MD Medical device compliant with EU Regulation no. 2017/745

TRIAL FRAME OCULUS UB6

7.6

For optometrical refractions - Tabo System

- 1 - PD adjustment, right and left eye
- 2 - Bridge height adjustment
- 3 - Corneal vertex distance adjustment
- 4 - PD measurement scale
- 5 - Scale for reading the axis position and base position
- 6 - Scale for measuring the corneal vertex distance (VD)
- 7 - Temple inclination adjustment
- 8 - Adjustable temple tips
- 9 - Temple length slide adjustment
- 10 - Knob for setting the axis position and base position
- 11 - Lens holder
- 12 - Nose pad on rack



NEW

Ref. 08238 Kr 12 243,-

Lightweight, non-allergic structure, suitable for all standard Ø 38 mm trial lenses.


- Determination of the astigmatism axis according to the TABO system
- Back vertex distance measurement gauge
- Independent monocular adjustments for optical centres
- Vertical and projection bridge adjustments
- Temple length and inclination adjustable
- PD adjustable range 48 to 80 mm
- Right and left PD adjustable independently from 24 to 40 mm
- Bridge height adjustable by up to 15 mm
- Temple adjustable in length from 98 to 135 mm and in inclination from +6° to -15°
- **6 lens holders each eye: 4 external rotating by 360° and 2 fixed internal**
- Axis rotation is continuously adjustable, in increments of 5°.

Weight: 63 g (without polarizing filters)

MD MD Medical device compliant with EU Regulation no. 2017/745

7.6 MENISCUS TRIAL LENS SET

236 meniscus \varnothing 38 mm trial lenses in carrying case



Ref. 08251 **Kr 23 250,-**

236 meniscus lenses in carrying case.

Lenses in optical glass with refraction index of 1.523.

- 35 pairs of positive sphere lenses from +0.12D to +20.00D
- 35 pairs of negative sphere lenses from -0.12D to -20.00D
- 19 pairs of positive cylinder lenses from +0.12D to +6.00D
- 19 pairs of negative cylinder lenses from -0.12D to -6.00D
- 10 prisms from 1 Δ to 10 Δ
- 10 supplementary lenses:

Maddox Rod white (2 pieces) - Grid - Pinhole - Occluder - Stenopeic slot - Red filter - Green filter - Clear glass - Opaque glass.

All the inside lens curves are the same to guarantee a consistent vertex distance.

The rims are made from polycarbonate and are coloured red for minus powers and black for plus powers and grey for prisms and supplementary lenses.

Case dimensions: 530 x 400 x 100 mm

Weight: 5.1 kg



MD MD Medical device compliant
with EU Regulation no. 2017/745

REGULAR TRIAL LENS SET

7.6

Trial lenses \varnothing 38 mm in carrying case



Ref. 08260 Kr 8755,-

232 plain lenses in carrying case with removable tray.
Lenses in optical glass with refraction index of 1.523.
Thin lenses with traditional flat bi-concave and bi-convex curves.

- 34 pairs of positive sphere lenses from +0.12D to +20.00D
- 34 pairs of negative sphere lenses from -0.12D to -20.00D
- 19 pairs of positive cylinder lenses from +0.12D to +6.00D
- 19 pairs of negative cylinder lenses from -0.12D to -6.00D
- 10 prisms from 1 Δ to 10 Δ
- 10 supplementary lenses:

Maddox Rod white - Pinhole (2 pcs.) - Occluder (2 pcs.) -
Stenopeic slot - Red filter - Green filter - Clear glass -
Opaque glass.

The rims are made from ABS plastic and are coloured red for minus powers and black for plus powers and grey for prisms and supplementary lenses.

Dimensions:

Case: 555 x 370 x 95 mm

Tray: 520 x 340 x 35 mm

Weight: 4.6 kg

NEW

ADDITIONAL LENSES

(not included in Ref. 08260)

Ref. 08260/0.25PRIS 0.25 Δ - 1 pc.

Ref. 08260/0.50PRIS 0.50 Δ - 1 pc.

Ref. 08260/0.75PRIS 0.75 Δ - 1 pc.

Ref. 08260/POL Linear polarized lens
45°-135° - 1 pair



MD MD Medical device compliant
with EU Regulation no. 2017/745

7.6 REGULAR TRIAL LENS SET

232 metal rim plain lenses Ø 38 mm in carrying case



Ref. 08261 **Kr 11 491,-**

232 plain lenses Ø 38 mm in optical glass with refraction index of 1.523.

Removable tray.

High dioptrical power thin lenses.

- 34 pairs of positive sphere lenses from +0.12D to +20.00D
- 34 pairs of negative sphere lenses from -0.12D to -20.00D
- 19 pairs of positive cylinder lenses from +0.12D to +6.00D
- 19 pairs of negative cylinder lenses from -0.12D to -6.00D
- 10 prisms from 1Δ to 10Δ
- 10 supplementary lenses:
 - Maddox Rod white - Pinhole (2 pcs.) - Occluder (2 pcs.) -
 - Stenopeic slot - Red filter - Green filter - Clear glass -
 - Opaque glass.

The rims are made of metal and are coloured red for minus powers and black for plus powers and grey for prisms and supplementary lenses.

Dimensions:

Case: 555 x 370 x 95 mm

Tray: 520 x 340 x 35 mm

Weight: 4.6 kg



MD MD Medical device compliant
with EU Regulation no. 2017/745

REGULAR TRIAL LENS SET

7.6

90 metal rim plain lenses Ø 38 mm in carrying case



Ref. 08262 Kr 5745,-

90 plain lenses Ø 38 mm in optical glass with refraction index of 1.523.

Removable tray.

- 17 pairs of positive sphere lenses from +0.25D to +10.00D
- 17 pairs of negative sphere lenses from -0.25D to -10.00D
- 8 pairs of negative cylinder lenses from -0.25D to -3.00D
- 7 supplementary lenses:

Maddox Red - Pinhole - Occluder - Stenopeic slot - Red filter - Green filter - Cross cylinder -0,25D.

The rims are made of silvered metal.

Dimensions:

Case: 370 x 280 x 90 mm

Tray: 345 x 250 x 30,5 mm

Weight: 2,8 kg




MD Medical device compliant with EU Regulation no. 2017/745



Ref. 08269 Kr 13 721,-

Suitcase for eye check/home care

Tablet not included



7.6 PROGRESSIVE TRIAL LENS SET

14 metal rim progressive AR coated lenses Ø 38 mm in carrying case



For best results, it is recommended to combine it with a trial lens set with lenses with a metal rim to ensure a larger field of vision

14 PROGRESSIVE LENSES

Ref. 08266 Kr 6019,-

Trial lens set with 14 progressive lenses Ø 38 mm in metal rims.

Used by the optician to let the customer experience progressive lenses before purchasing ophthalmic lenses.

AR coated lenses in organic material with refraction index of 1.56.

Short progressive corridor length of 12 mm

Addition power ranges included:

+1.00D, +1.25D, +1.50D, +1.75D, +2.00D, +2.50D, +3.00D.

For each addition power range there is one right lens and one left lens.

Laser engravings identify:

- Nasal position "N"
- Horizontal axis "-"
- Near fitting circle "O"

The rims are made of metal.

Removable tray.

Dimensions:

Case: 175 x 155 x 95 mm

Tray: 150 x 130 x 40 mm

Weight: 0.9 kg



MD Medical device compliant
with EU Regulation no. 2017/745

PRISM TRIAL LENS SET

7.6

76 prism Ø 38 mm trial lenses with AR coating in carrying case



76 PRISM LENSES

Ref. 08264 Kr 10 807,-

76 prism lenses Ø 38 mm in carrying case with removable wooden tray. Suitable for the measuring and correction methodology after H.-J. Haase (MKH).

AR coated lenses in optical glass with refraction index of 1.523.

Prism diopter range:

0.25Δ-0.50Δ-0.75Δ-1.00Δ-1.25Δ-1.50Δ-1.75Δ-2.00Δ-2.50Δ-3.00Δ-3.50Δ-4.00Δ-4.50Δ-5.00Δ-6.00Δ-7.00Δ-8.00Δ-9.00Δ-10.00Δ.

This set contains 4 lenses for each prism diopter, with 4 different prism inclinations.

A reference mark indicates the base (notch) and the vertex (dot) of the prism.

Thin lenses with high prism powers.

The rims are made from grey coloured ABS plastic.

Dimensions:

Case: 300x190x93 mm

Tray: 270x155x35 mm

Weight: 1.8 kg

MD MD Medical device compliant with EU Regulation no. 2017/745

7.6 TRIAL CLIPS - TRIAL LENS HOLDER

Trial Clip

Ref. 08229 **Kr 1422,-**

Monocular trial clip which may be used with almost all types of spectacles frames.

2 trial clips with 3 lens holders for \varnothing 38 mm trial lenses.

Ideal for low vision testing.

Weight: 9 g each





Trial Lens Holder for 4 lenses

Ref. 08344 **Kr 246,-**

Trial Lens Holder for 4 lenses

With 0°-45°-90°-135° and 180° axis marks

For \varnothing 38 mm trial lenses.



MD MD Medical device compliant
with EU Regulation no. 2017/745

DIGITAL PD METER

7.7

With PD adjustment slide



- 1 - PD/LD switch
- 2 - PD adjustment slide (left)
- 3 - Focusing distance selection lever
- 4 - Practitioner's viewpiece
- 5 - Viewpiece add-power switch
- 6 - PD adjustment slide (right)
- 7 - ON/OFF switch
- 8 - Right eye measurement display
- 9 - PD distance
- 10 - Left eye measurement display
- 11 - Eye shield lever (mono/binocular switch)

Ref. 04966 **Kr 5745,-**


With PD adjustment slides.
Additional lens +2.00D switchable on the PD measurement ocular.

- Internal light with high brilliancy LED
- Measurement range: 45-82 mm (0,5 mm increments)
- Measurement: monocular and binocular
- Viewing distance: from 30 cm to infinity
- Power off: automatically after 1 minute without use

Power supply: 2 batteries (1.5V AA) not included

Dimensions: 221x166x63 mm

Weight: 0,66 kg (with batteries)



MD MD Medical device compliant with EU Regulation no. 2017/745

7.7 REFRACTOR/PHOROPTER - VISION TESTER

Ophthalmic instrument used to investigate sight problems and binocular function



Ref. 08249 Kr 35 978,-

The vision tester is an instrument that is used to examine refractive problems. It contains most of the lenses used in the trial lens set including the Resley prisms. Pupillary distance adjustment: from 48 to 75 mm in increments of 1 mm.

- Antireflection lenses
- Sphere lenses: from +16.75D to -19.00D (in increments of 0.25D)
- Cylindrical lenses: from -0.00D to -6.00D in increments of 0.25D (extendible to -8.00D with additional -2.00D auxiliary lens)
- Cylindrical axis: from 0 to 180° (in increments of 5°)
- Cross cylinders: ±0.25 D
- Rotating prisms: from 0Δ to 20Δ in increments of 1Δ

Dimensions: 345x350x155 mm
Weight: approx. 4.5 Kg

- Supplementary lenses:
 - (O) Open lenses
 - (R) Retinoscopy lenses +1.50 D
 - (P) Polarized lenses: (45° left eye/135° right eye)
 - (WMV) or (RMV) Maddox, vertical: left eye - white; right eye - red
 - (WMH) or (RMH) Maddox, horizontal: left eye - white; right eye - red
 - (GL) Green lens
 - (RL) Red lens
 - (+.12) Sphere lens + 0.12 D
 - (PH) Pinhole
 - (10 Δ I) o (6Δ U) Dissociated prisms: left eye 10 (Δ) I - right eye 6 (Δ) U
 - (±.50) Fixed cross cylinders
 - (OC) Occluder

MD MD Medical device compliant with EU Regulation no. 2017/745

SELF BALANCING ARM

7.7

For Phoropter and Projector



ACCESSORIES INCLUDED WITH REF. 08249:

- Near vision test card with rod
- Protective cover
- Face shield with protective lens.

Ref. 08243


Kr 10 900,-

With wall-mounted support for vision tester, complete with chart projector mounting.



7.7 MULTIFUNCTION RULER

9 measurement functions in one tool



NEW

Ref. 04961 Kr 437,-

This multifunction ruler was designed to include 9 of the most useful functions for opticians for the correct fitting of lenses to the frame. Light, compact and easily transportable, it can be used both in the practice and during home visits.

1 pc.

- 1 Ruler: graduated scale from 0 to 200 mm.
- 2 Wrapping angle: curvature of the front of the frame is designed to keep the lenses as much as possible equidistant for all viewing positions on the horizontal plane. Most optical frames have a curve of about 5.
- 3 Lens diameter: Minimum diameter that the lens must have before being cut.
- 4 Progression channel height: millimeter reference for measuring the height of the canal progression.

- 5 Pantoscopic angle: angle of inclination of the front with respect to the plane determined by the temples. The front must be tilted down towards the face of the patient in order to ensure as much as possible the equidistance of the lenses to the eye for all positions of view. Most optical frames have a pantoscopic angle of about 7.
- 6 Interpupillary distance: Graduated scale for measurement of the two half distances or of the total pupillary distance.
- 7 Height of the bifocal lens progression channel: Height of the progression channel that joins the distance optical zone with the reading area zone.
- 8 Corneal apex distance - lens: distance from internal surface of the ophthalmic lens at the corneal apex, usually between 12 and 14 mm.
- 9 Base of the frame: Double measurement scale to check the base of different frames size and caliber.


TEST ROOM ACCESSORIES

7.7

PD Scales

Ref. 04952 **Kr 54,-**
 Plastic PD ruler to check children's PD. Measuring range: 100 mm


 1 pc.



Transparent section with millimeter scale for fast and quick check. Highly flexible, rounded edges.


Ref. 04959 **Kr 68,-**
 PLASTIC PD SCALE WITH EYE OCCLUDER
 Measuring range: 150 mm - 1 pc.

 1 pc.



Ref. 04954 **Kr 52,-**
 PLASTIC - 170 MM
 Highly flexible, scale non-erasable - 1 pc.

 1 pc.



Adhesive film for bifocal and progressive heights



FITTING HEIGHT MEASURING TAPE

Ref. 04963 **Kr 1436,-**

Static cling film that quickly and easily determines bifocal seg height and progressive fitting cross placement on any frame and the depth of the groove.

- 1 mm increments.

Roll 500 pcs.



7.7 READING CHARTS FOR 3 METRES DISTANCE TESTS

Flexible plastic material



Ref. 0842 Kr 957,-

Reading Chart made of washable flexible plastic material.

- Testing distance approx. 3 metres
- Can be fixed to the wall or used with a light box for rear illumination
- SLOAN letters
- LogMAR notation
- Line progression in LogMAR, decimal, metres, foot.


Dimensions: 230x360 mm
Weight: 65 g

Ref. 08426 Kr 820,-

Reading Chart made of washable flexible plastic material.

- Testing distance approx. 3 metres
- Can be fixed to the wall or used with a light box for rear illumination
- LEA symbols, ideal for children from the age of 3 years
- LogMAR notation
- Line progression in LogMAR, decimal, metres, foot.


Dimensions: 230x360 mm
Weight: 65 g



MD MD Medical device compliant with EU Regulation no. 2017/745

READING CHARTS FOR 3 METRES DISTANCE TESTS **7.7**

Flexible plastic material



Ref. 08427 Kr 820,-

Reading Chart made of washable flexible plastic material.

- Testing distance approx. 3 metres
- Can be fixed to the wall or used with a light box for rear illumination.
- "E" letters specific for examinations on children from 4 to 6 years old.
- Lines progression: 10 ft. and 20 ft.

Dimensions: 230x360 mm

Weight: 65 g

NEW

Ref. 08424 Kr 280,-

Letter E in paper version.

Pack of 30 sheets size 210x297 mm in black

cardboard with 2 letter E

Letter size: 90x160 mm



MD MD Medical device compliant
with EU Regulation no. 2017/745

7.7 OPTOTYPES FOR NEAR VISION

For reading test



OPTOTYPE "ACTIVE"
Ref. 08421GB Kr 109,-
 Paper optotype for reading test for children.
 Dimensions: 150 x 210 mm.
 2 pcs.



OPTOTYPE "CENTRO STYLE"
Ref. 08422GB Kr 109,-
 Paper optotype for reading test for adults.
 Dimensions: 165 x 290 mm.
 2 pcs..



OPTOTYPE "AIRPORT"
Ref. 08423GB Kr 109,-
 Paper optotype for reading test for adults.
 Dimensions: 150 x 210 mm. 2 pcs.

Schirmer Test - Fluorescein Sodium - Easy Mark



SCHIRMER TEST

Ref. 08781 Kr 478,-
 Graduated strips for dry eye test. Measurement of tear volume and production. 100 pcs. (strips in single sterile pack)

FLUORESCEIN SODIUM

Ref. 08783 Kr 383,-
 Evaluation of integrity of the ocular surface, G.P. hard contact lens fitting and applanation tonometry. Each strip is impregnated with 1.0 mg of fluorescein sodium.

100 pcs.
 (strips in single sterile pack)



EASY MARK

Ref. 08793 Kr 275,-
 Sterile solution based on Riboflavin for the diagnosis of the ocular surface and staining of the tear film. Certified for applanation tonometry and the application of hard and soft contact lenses. 10 single-dose vials - 0,5 ml each



7.8 ACCESSORIES FOR CONTACT LENSES

Cases For soft and G.P. hard lenses



Suitable also for:

- lenses used in orthokeratology.
- lenses with a diameter between 8.8 and 11 mm

Ref. 08747
Kr 168,-
6 pcs. assorted



Ref. 08730
Kr 205,-
6 pcs. assorted



MD MD Medical device compliant
with EU Regulation no. 2017/745

Suckers Holders and Tweezers



SUCKER HOLDERS
 Indicate for semi-rigid contact lenses

- Ref. 08764 Kr 246,-
- Ref. 08765 Kr 246,-

12 pcs./ref.



WITH HOLLOW CENTER
 Indicate for soft contact lenses

- Ref. 08767 Kr 273,-
- Ref. 08766 Kr 273,-

12 pcs./ref.

TWEEZERS
 Ref. 08772 Kr 191,-
 With case
 12 pcs. assorted

TWEEZER
 Ref. 08770 Kr 177,-
 With case
 6 pcs.



7.8 DISPOSABLE TISSUES





VIDEO



Protection during the visit




Ref. 06983 Kr 207,-
Disposable paper tissues for front rests and for optical equipment.
Box of 200 pcs.



Ref. 06984 Kr 191,-
Disposable paper tissues for front rests for phoropter and optical equipment. Box of 200 pcs.




Ref. 06980 Kr 118,-
Disposable paper tissues for chin rests and for optical equipment. Box of 500 pcs.



LAMINATED POSTER "THE EYE"

7.8

Anatomical Image of the eye



NEW

Ref. 11105 Kr 266,-

Poster illustrates the different parts of the eye i.e muscle, crystalline, optical nerve, lacrimal gland, net, cornea, eyelashes corps.

Dimensions: 57,5x42,5 cm

1 pc.



Kontaktinformasjon:
Rodenstock Utstyersavdelingen
Tlf: 915 95 231
E-post: utstyersavdelingen@rodenstock.no

Alle priser er eksklusiv merverdiavgift,
og er gyldige til 31.12.2025.

Vi tar forbehold om trykkfeil og eventuelle endringer.